



## BROMO THYMOL BLUE PAPER (10 BKS)

CAT. NO. : PH0011  
GRADE : LR  
CAS NO. : \_  
FORMULA : \_  
M.W. : \_

Orange-yellow strips

Sensitivity

To pass the test.

### pH Indicator Papers – Important Usage Information

Slight discoloration of the dry, unused indicator papers may occur due to their high sensitivity to environmental conditions and pH exposure. This does **not** affect the performance or reliability of the papers in pH determination.

When testing **very weakly buffered or unbuffered solutions**, an *indicator error* may occur. This can be minimized by following the method below:

1. Press the indicator strip against the inner wall of a clean test tube.
2. Fill the test tube with the test solution until the strip is fully immersed up to its top edge.
3. After **30 seconds to 1 minute**, compare the color of the strip with the pH scale **through the wall of the test tube**.

DATE: 01.07.2025

**Remarks:** The above-mentioned batch material is compliance with the standard specifications of High Purity Laboratory Chemicals Pvt. Ltd.

