



+91.22.22057323
+91.84879 47881

Product Specification

UNIVERSAL INDICATOR PAPER PH RANGE 1-14

CAT. NO. : PH0097
GRADE : LR
CAS NO. : _
FORMULA : _
M.W. : _

ISO 9001:2015 CERTIFIED

Sensitivity

Passes test



♦ pH-Indicator Paper – Usage Note

- Slight discoloration of dry, unused papers may occur due to their high sensitivity to pH and ambient exposure. This does **not** affect their performance in pH determination.
- In weakly buffered or unbuffered solutions, **indicator error** may arise. To minimize this:
 - Adhere the strip to the **inner wall of a test tube**.
 - Fill the tube up to the top edge of the paper with the test solution.
 - Wait **30 seconds to 1 minute**.
 - Compare the paper's color with the reference scale **through the glass**.

DATE: 01.07.2025

Remarks: The above-mentioned batch material is compliance with the standard specifications of High Purity Laboratory Chemicals Pvt. Ltd.



Customercare@hplcindia.com



hplcindia.com